MARSDEN

MG-4800 Hand Grip Dynamometer



- 20kg/40kg/80kg spring resistance
- · Ideal for occupational therpists

- · Max mode can measure maximum grip force
- AVG mode can give the average force

Marsden MG-4800:

Hand Grip Dynamometer

- 20kg/40kg/80kg spring resistance
- Ideal for occupational therpists
- Max mode can measure maximum grip force
- AVG mode can give the average force

- Can be used in the home environment
- Easy to operate
- Enables users to develop strength in their hand and forearms
- Covered by Marsden 4 year warranty



Free standing design

The MG-4800 is designed to help rehabilitation of people with disabilities or injuries by enabling them to develop strength of their hands and forearms.

We have found the MG-4800 to be incredibly popular with occupational therapists, rehabilitation doctors and healthcare assessment clinics. Also, for continual home use by patients.

Help to build muscle

Regular use of the MG-4800 will assist with a gradual build of the muscles and grip, providing the user with the opportunity to regain critical abilities.

The hand grip dynometer is sold as standard with three different kinds of resistance springs at 20kg, 40kg and 80kg. These can be changed as and when the muscle improves and milestones are achieved.

Easy to use

The MG-4800 is a very single and easy to use device. Ideal for use by medical professionals and people at home wanting to improve their strength.

With three simple buttons, the device is easy to operate and it can easily switch between modes. The resistance springs are easy to chance and can be done at home without difficulty.

Resistance:	20kg/40kg/80kg	Modes:	Max mode, AVG
Warranty:	4 Years	Weight::	0.6kg