

MARSDEN

T-120 Non-Contact Forehead and Ear Thermometer



- MDD Approved
- 3-in-1 tympanic and non-contact thermometer
- Simple and easy to use
- Measures temperature in °C or °F
- 20 measurement memory
- 1 year warranty

Marsden T-120:

Tympanic and Forehead Thermometer

- Accurate to 0.1°C
- MDD Approved
- Forehead, ear and object infrared thermometer
- Clear and easy to read digital reading
- Simple point and press for quick thermometer reading
- Measuring distance between 1 to 3cm
- Measures temperature in °C or °F
- Auto-power off to prolong battery life
- Replaceable batteries to aid portability
- Can be used to measure bath water, baby milk and surfaces
- Detects signs of virus in seconds
- An alert will signal when higher than 37.8 °C



What is the T-120?

The T-120 is an MDD Approved digital infrared thermometer that makes taking a persons temperature incredibly simple and easy to use.

This is a handheld device which is ultra-lightweight and easily portable. Temperature readings can be provided in either °C or °F and it is accurate to 0.1°C to ensure that people are monitored closely and safely.

How to use the T-120?

The T-120 can be used either by pointing at an individuals forehead, alternatively remove the cap, place on a cover and use as a tympanic thermometer.

Once a temperature has been registered on the device, it will show nice and clear on the display. If an individuals temperature is above 37.8 °C, the display will turn red and an alert will beep.

Where can it be used?

Place the thermometer 1 to 3cm away from the forehead and press the button, or apply a cover and place in their ear for a high degree of accuracy.

The secondary function of the T-120 is to monitor the temperature of surfaces such as bath water, baby milk and alternative surfaces. Easy to toggle between by a simple button.

Memory:	20 items	Power Supply:	Replaceable Batteries (2 x AAA batteries)
Graduations:	0.1°C	Battery Life:	Upto 300 measurements
Functions:	Unit Switchable, Portable	Approvals:	MDD Approved
Dimensions:	146mm x 37mm x 39mm	Optional Extras:	Covers
Weight of Scale:	0.68kg	Warranty Period:	1 Year