

V-100 Best-Selling Veterinary Scale



- Capacity: 100kg
- Accurate from 50g
- Portable and lightweight pet scale

- Simple and easy to use
- Tare to remove unwanted weight
- Covered by 1 year warranty

Marsden V-100:

Best-Selling Veterinary Scale

- Capacity to 100kg
- Accurate from 50g
- Portable scale for weighing pets
- Lightweight and easy to transport
- Automatic hold for accuracy
- Perfect size for cats and dogs

- Weigh a wide range of animals with ease
- Wipe down top
- Easy to clean and keep hygienic
- Simple and easy to use
- Tare to remove unwanted weight
- Covered by 1 year warranty



What is the V-100?

The V-100 is Marsden's best-selling veterinary scale. It is a 100kg capacity scale, perfect for weighing small to medium sized pets, particularly cats and dogs.

This scale features from Tare and Hold functions and a clear digital weight display to ensure that all users, whether you be in a veterinary surgery or at home, are able to obtain accurate weights.

Easy to clean

In order to help make the V-100 incredibly diverse and suitable for veterinary surgeries, we have designed the product to be wiped down.

Even through the device does not hold an IP rating, the platter can be wiped down without causing damage to the scale. Making it a clean and hygienic addition to any surgery.

Lightweight and portable

The V-100 has been designed with the ease of the user in mind. Whether it be within a veterinary surgery or at home, this is the ideal scale.

Weighing only 8kg, the scale is easy to transport around the veterinary surgery, from room to room. Or when used at home it is slimline enough to be stored in small and compact spaces.

Capacity:	100kg	Power Supply:	Mains or Replaceable (2 x AAA batteries)
Graduations:	50g	Battery Life:	40 hours of continuous use
Functions:	Hold, Portable, Tare	Approvals:	Non-Approved
Dimensions:	715mm x 410mm x 40mm	Optional Extras:	N/A
Weight of Scale:	8kg	Warranty Period:	1 Year