

MARSDEN

V-110 Large Veterinary Scale



- Capacity: 150kg
- Accurate from: 100g
- Portable and lightweight pet scale
- Simple and easy to use
- Tare to remove unwanted weight
- Covered by 1 year warranty

Marsden V-110:

Large Veterinary Scale

- Capacity to 150kg
- Accurate from 100g
- Portable scale for weighing pets
- Lightweight and easy to transport
- Automatic hold for accuracy
- Perfect size for cats and dogs
- Weigh a wide range of animals with ease
- Wipe down top
- Easy to clean and keep hygienic
- Simple and easy to use
- Tare to remove unwanted weight
- Covered by 1 year warranty



What is the V-100?

The V-110 is Marsden's most affordable high capacity veterinary scale. This portable, one-piece animal scale is suitable for weighing medium to larger animals.

The capacity of 150kg ensures that the scale will not be overloaded easily by your pets, making this a long-lasting and reliable weighing scale. Ideal for use in vet surgeries or at home.

Easy to clean

In order to help make the V-110 incredibly diverse and suitable for veterinary surgeries, we have designed the product to be wiped down.

Even though the device does not hold an IP rating, the platter can be wiped down without causing damage to the scale. Making it a clean and hygienic addition to any surgery.

Lightweight and portable

The V-110 has been designed with the ease of the user in mind. Whether it be within a veterinary surgery or at home, this is the ideal scale.

Weighing only 5.5kg, the scale is easy to transport around the veterinary surgery, from room to room. Or when used at home it is slimline enough to be stored in small and compact spaces.

Capacity:	150kg	Power Supply:	Mains or Replaceable
Graduations:	100g	Battery Life:	40 hours of continuous use
Functions:	Hold, Portable, Tare	Approvals:	Non-Approved
Dimensions:	900mm x 500mm x 62mm	Optional Extras:	N/A
Weight of Scale:	5.5kg	Warranty Period:	1 Year