MARSDEN

V-180 Large Veterinary Scale



Capacity: 180kg

Accurate from: 100g

Lightweight base and stainless steel platter

- · Simple and easy to use
- · Tare to remove unwanted weight
- Covered by 1 year warranty

Marsden V-180:

Large Veterinary Scale

- Capacity to 180kg
- Accurate from 100g
- Lightweight and easy to transport
- Automatic hold for accuracy
- Perfect size for cats and dogs
- Weigh a wide range of animals with ease

- Wipe down top stainless steel platter
- Easy to clean and keep hygienic
- Simple and easy to use
- Tare to remove unwanted weight
- Powered by 4 x AA batteries
- Covered by 1 year warranty



What is the V-180?

The V-180 is an affordable veterinary scale for weighing larger pets up to 180kg in weight. This scale is large, lightweight and has a stainless steel platter.

The scale is simple to use and features a large and bright LCD display to ensure that all users can easily see the weight reading, regardless of which environment they are using the scale in.

Easy to clean

To ensure that longevity of this animal weighing scale, we designed this model to be have a stainless steel construction, meaning it can be cleaned easily.

This means that animals can stand on the scale when they are wet and muddy, with no detrimental effect to the weighing scale and it can be wiped down by the veterinary staff and it is good to go.

Highly portable

We understand that whether the scale is been used at home or within a veterinary surgery, it is highly likely that the scale will need to be transported and moved from time to time.

For this reason, the V-180 has been designed to be super lightweight at only 6.2.kg and has a carry handle to make transportation easier.

Capacity:	180kg	Power Supply:	Mains or Replaceable (4 x AA batteries)
Graduations:	100g	Battery Life:	40 hours of continuous use
Functions:	Hold, Portable, Tare	Approvals:	Non-Approved
Dimensions:	900mm x 500mm x 25mm	Optional Extras:	N/A
Weight of Scale:	6.2kg	Warranty Period:	1 Year