MARSDEN

V-22 Veterinary Scale



- Capacity: 20kg
- Accurate from 5g
- Portable vet scale with removable tray
- · Simple and easy to use
- Tare to remove unwanted weight
- Covered by 1 year warranty

Marsden V-22:

Small Veterinary Scale

- Capacity to 20kg
- Accurate from 5g
- Portable scale for weighing small pets
- Lightweight and easy to transport
- Automatic hold for accuracy
- Perfect size for small cats and dogs

- Weighs small animals such as rabbits and hamsters with ease
- Wipe down top
- Easy to clean and keep hygienic
- Simple and easy to use
- Tare to remove unwanted weight
- Covered by 1 year warranty



What is the V-22?

The V-22 is a premium veterinary scale with a removable tray which can be used as a stretcher. It's suitable for weighing domestic pets in vet surgeries and at home.

It can accommodate small animals, such as cats and dogs, providing weight readings to an accuracy of 5g and with a capacity of 20kg. Highly popular due to the small and compact size of the scale.

Easy to clean

In order to help make the V-22 incredibly diverse and suitable for veterinary surgeries, we have designed the product to be wiped down.

Even through the device does not hold an IP rating, the platter can be wiped down without causing damage to the scale. Making it a clean and hygienic addition to any surgery.

Lightweight and portable

The V-22 has been designed with the ease of the user in mind. Whether it be within a veterinary surgery or at home, this is the ideal scale.

Weighing only 2.5kg, the scale is easy to transport around the veterinary surgery, from room to room. Or when used at home it is slimline enough to be stored in small and compact spaces.

Capacity:	20kg	Power Supply:	Mains or Replaceable
Graduations:	5g	Battery Life:	40 hours of continuous use
Functions:	Hold, Portable, Tare	Approvals:	Non-Approved
Dimensions:	585mm x 335mm x 115mm	Optional Extras:	N/A
Weight of Scale:	2.7kg	Warranty Period:	1 Year