

MARSDEN

V-24 Veterinary Scale



- Capacity: 20kg
- Accurate from 10g
- Portable vet scale with removable tray
- Simple and easy to use
- Tare to remove unwanted weight
- Covered by 1 year warranty

Marsden V-24:

Portable Veterinary Scale

- Capacity to 20kg
- Accurate from 10g
- Portable scale for weighing small pets
- Lightweight and easy to transport
- Automatic hold for accuracy
- Perfect size for small cats and dogs
- Weighs small animals such as rabbits and hamsters with ease
- Wipe down top
- Easy to clean and keep hygienic
- Simple and easy to use
- Tare to remove unwanted weight
- Covered by 1 year warranty



What is the V-24?

The V-24 is an accurate, easy to use portable veterinary scale. The modern, professional scale provides digital weight readings to the nearest 10g.

This is the perfect addition to any veterinary surgery or for individuals that are looking to closely monitor the weight of their pet at home. Highly robust device with a built in overload protection.

Easy to clean

In order to help make the V-24 incredibly diverse and suitable for veterinary surgeries, we have designed the product to be wiped down.

Even though the device does not hold an IP rating, the platter can be wiped down without causing damage to the scale. Making it a clean and hygienic addition to any surgery.

Lightweight and portable

The V-24 has been designed with the ease of the user in mind. Whether it be within a veterinary surgery or at home, this is the ideal scale.

Weighing only 2.5kg, the scale is easy to transport around the veterinary surgery, from room to room. Or when used at home it is slimline enough to be stored in small and compact spaces.

Capacity:	20kg	Power Supply:	Mains or Replaceable (3 x AAA Batteries)
Graduations:	10g	Battery Life:	60 hours of continuous use
Functions:	Hold, Portable, Tare	Approvals:	Non-Approved
Dimensions:	609mm x 326mm x 91mm	Optional Extras:	N/A
Weight of Scale:	1.55kg	Warranty Period:	1 Year