MARSDEN

V-250 Large Veterinary Scale



· Capacity: 250kg

Accurate from: 100g

Lightweight base and stainless steel platter

- · Simple and easy to use
- · Tare to remove unwanted weight
- Covered by 1 year warranty

Marsden V-250:

Large Veterinary Scale

- Capacity to 250kg
- Accurate from 50g
- Automatic hold for accuracy
- Perfect size for cats and dogs
- Weigh a wide range of animals with ease
- Wipe down top stainless steel platter

- Easy to clean and keep hygienic
- Simple and easy to use
- Tare to remove unwanted weight
- Optional rubber matt available
- Powered by 6 x AA batteries
- Covered by 1 year warranty



What is the V-250?

The V-250 is a large veterinary scale with a capacity of 250kg. This scale is has a supersized stainless steel base which is suitable for weighing even the largest of animals.

An optional rubber mat is available, which simply lays over the top of the stainless steel platter, to ensure that animals are safe and secure whilst being weighed on the V-250.

Easy to clean

To ensure that longevity of this animal weighing scale, we designed this model to be have a stainless steel construction, meaning it can be cleaned easily.

This means that animals can stand on the scale when they are wet and muddy, with no detrimental effect to the weighing scale and it can be wiped down by the veterinary staff and it is good to go.

Heavy duty

We understand here at Marsden that you when you are weighing large animals, you want the equipment to not only be accurate, but be long-lasting and reliable.

The V-250 has been constructed to the highest quality and that ensures that the device will be long-lasting and reliable for many years to come, regardless of if it is been used at home or the surgery.

Capacity:	250kg	Power Supply:	Mains or Replaceable (6 x AA batteries)
Graduations:	100g	Battery Life:	40 hours of continuous use
Functions:	Hold, Portable, Tare	Approvals:	Non-Approved
Dimensions:	1090mm x 505mm x 70mm	Optional Extras:	Rubber Mat
Weight of Scale:	18kg	Warranty Period:	1 Year