

MARSDEN

VT-250 Veterinary Table Scale



- Capacity: 250kg
- Accurate from: 100g
- Vet table designed for the consulting room
- Simple and easy to use
- Tare to remove unwanted weight
- Covered by 1 year warranty

Marsden VT-250:

Large Veterinary Consulting Room Table Scale

- Capacity to 250kg
- Accurate from 50g
- Table top scale, designed for use in the consulting room
- Automatic hold for accuracy
- Perfect size for cats and dogs
- Weigh a wide range of animals with ease
- Wipe down top stainless steel platter
- Easy to clean and keep hygienic
- Simple and easy to use
- Tare to remove unwanted weight
- Sold as standard with a rubber mat
- Covered by 1 year warranty



What is the VT-250?

The VT-250 is a large veterinary scale with a capacity of 250kg. This scale is built within a table top, making it the ideal weighing solution for veterinary surgeries.

Sold as standard with a rubber mat. Simply lay the rubber mat over the top of the stainless steel platter, to ensure that animals are safe and secure whilst being weighed on the VT-250.

Easy to clean

To ensure that longevity of this animal weighing scale, we designed this model to be have a stainless steel construction, meaning it can be cleaned easily.

This means that animals can stand on the scale when they are wet and muddy, with no detrimental effect to the weighing scale and it can be wiped down by the veterinary staff and it is good to go.

Heavy duty

We understand here at Marsden that you when you are weighing large animals, you want the equipment to not only be accurate, but be long-lasting and reliable.

The VT-250 has been constructed to the highest quality and that ensures that the device will be long-lasting and reliable for many years to come, regardless of if it is been used at home or the surgery.

Capacity:	250kg	Power Supply:	Mains or Replaceable (6 x AA batteries)
Graduations:	100g	Battery Life:	40 hours of continuous use
Functions:	Hold, Portable, Tare	Approvals:	Non-Approved
Dimensions:	1160mm x 550mm x 830mm	Supplied With:	Rubber Mat
Weight of Scale:	35kg	Warranty Period:	1 Year