

# MARSDEN

## B-200 Small Veterinary Scale



- Various capacities available
- Accurate from 5g
- Portable scale for weighing small pets
- Separate indicator with a wall bracket
- Tare to remove unwanted weight
- Covered by 1 year warranty

# Marsden B-200:

## Small Veterinary Scale

- Various capacities available
- Accurate from 5g
- Portable scale for weighing small pets
- Separate indicator with a wall bracket
- Our most popular bench top scale for small animals
- Perfect size for small cats and dogs
- Weighs small animals such as rabbits and hamsters with ease
- Wipe down stainless steel platter
- Easy to clean and keep hygienic
- Simple and easy to use
- Tare to remove unwanted weight
- Covered by 1 year warranty



### What is the B-200?

The B-200 is an ultra-lightweight table top weighing scale, suitable for weighing small pets such as hamsters, rabbits, cats and small dogs.

The scale features from a small stainless steel platter and an indicator on a flying lead which can be mounted on a wall if required. Very simple and easy to use with only five function keys.

### Easy to clean

In order to help make the B-200 incredibly diverse and suitable for veterinary surgeries, we have designed the product to feature from a stainless steel top plate.

Even though the device does not hold an IP rating, the stainless steel platter can be wiped down without causing damage to the scale. Making it a clean and hygienic addition to any surgery.

### Customise to your needs

The stand out thing about the B-200 is that it is a scale with many options available, so this can be tailored to your requirements, whatever they may be.

Various capacities are available, ranging from 15kg to 60kg. Depending on the capacity you choose, the accuracy will adjust, making this a customisable device, to suit your needs.

<b>Capacity:</b>	15kg/30kg/60kg	<b>Power Supply:</b>	Mains or Replaceable (6 x AA batteries)
<b>Graduations:</b>	5g/10g/20g	<b>Battery Life:</b>	100 hours of continuous use
<b>Functions:</b>	Hold, Portable, Tare	<b>Approvals:</b>	Non-Approved
<b>Dimensions:</b>	315mm x 305mm x 55mm	<b>Optional Extras:</b>	N/A
<b>Weight of Scale:</b>	4.1kg	<b>Warranty Period:</b>	1 Year